



# AMADEO

X-RAY SYSTEMS

## Portable High Frequency X-ray Source



### PSM-PD

#### Product Composition

PSM-PD consists of controller, MBK, beam defining clipper

#### Applications

PSM-PD is mainly applied to mobile DR, military DR and vehicle-carried DR. cle-carried DR.

#### Feature

- Wireless communication for remote monitor
  - Wide range AC input with APFC integrated (100~240VAC)
  - Soft switching technology, switching frequency is up to 250kHz
  - Intellectualized software system, automatic filament adjustment, dual-filament selection and self-detection capabilities
- Operation interface can be divided into human and veterinary, support for APR data configuration for different areas/parts

Product Model		PSM-PD5.6CPE
Product Name		Portable High Frequency High Voltage Generator
Input	Max Working Frequency	≥250kHz
	Phase/Voltage Range	100VAC~240VAC/48VDC
	Frequency Range	50Hz/60Hz±1Hz
Output	Nominal/Max Power	5.6kW/5.6kW
	kV Range	40kV~125kV
	kV Ripple	≤3%@100kV
	kV Rising Time	≤1ms@100kV(10%~90%)
	Accuracy of kV	±5% or ≤4kV
	mA Range	5mA~100mA
	Accuracy of Tube Current	±(5%+1mA)@≥10ms ±20%@ < 10ms
	mAs Range	0.1~320mAs@220VAC
		0.1~200mAs@110VAC
		0.1~100mAs@48VDC
	mAs Accuracy	±(5%+0.2mAs)@ > 5ms ±(10%+0.2mAs)@≤5ms
	Exposure Time Range	2ms~10000ms
	Accuracy of Exposure Time	±(2%+0.2ms)@ > 5ms ±(5%+1ms)@≤5ms
Tube Parameter	Tube Model	KL10-0.6/1.8-125
	Norminal Focal Spot Size	0.6mm/1.8mm
	Max Anode Heat Storage Capacity	20kJ
	Anode Disk Target Angle	15°
Communication		WIFI/RS232
Appearance	Dimensions (L*W*H)	383*240*215 mm
	Weight	15kg

## Dentem Sağlık Hizmetleri

Aydınevler Mah. Oğuzhan Cad. Ata Plaza  
No: 1-3 Kat:3 Kapı No: 31 Maltepe / **İstanbul**

+ 90 (216) 504 00 47  
info@dentem.com.tr  
www.dentem.com.tr



**ÜTS**  
ÜRÜN TAKİP SİSTEMİ  
ÜTS NO: 06973858210447 (GS1)

Suzhou Powersite Electric Co., Ltd.  
Floor 1st, 2nd, 4th, buildgin 5  
No. 188 Funchunjiang Road,  
Suzhou New District, Suzhou  
PEOPLE'S REPUBLIC OF CHINA